


Gap.Jumper:

a probabilistic approach for single nucleotide variant quality assessment from samples, replicates and software results

Pawel Rosikiewicz^{1*}, Frédéric G Masclaux^{1,2}, Tania Wyss¹,
Frédéric Schütz², Marco Pagni², Ian R. Sanders¹

¹ Department of Ecology and Evolution, University of Lausanne, Bâtiment Biophore, Lausanne, 1015, Switzerland
² Vital-IT, SIB, Swiss Institute of Bioinformatics, Bâtiment Génomode, Lausanne 1015, Switzerland,
 * mail: prosikie@unil.ch




Next generation sequencing (NGS) allows screening of genetic polymorphisms in samples with a high genetic polymorphisms such as samples of carcinoma or from organisms with different ploidy (e.g. pathogenic fungi)

Problems are:


- limited coverage
- missing data
- sequencing errors
- different results obtained with different software's
- technical and biological differences between replicates

Consequently, researchers are faced with data containing a large number of apparently variable positions that need to be confirmed with independent experimental approach


Gap.jumper allows integration of variant calling data obtained from different samples, replicates and software results




Gap.Jumper allows estimating uncertainty associated with each nucleotide based on available empirical data



Evidence scores can be used to remove potential errors or to rank positions based on their quality



Accuracy improves with increasing number of pooled replicates



Evidence scores can be used to estimate uncertainty of polymorphisms detected between different samples

Application example:


Identification of the fungal clonal lines

Goal:
To identify which fungal strains are clonal offspring produced from a common parent (benchmarking studies with known classes)

Methods:
We genotyped 42 fungal strains (RAD-seq) - three replicates of each strain were sequenced.

The replicates of each strain were used to build a consensus (DST-based approach), which was compared to consensus built for other strains

Two group of strains which were previously identified as clonal lines were used as control (black and red triangles and black and grey squares).



CONCLUSIONS:

- validates SNPs accurately in sets with a relatively small number of replicates (2-6)
- handles missing information easily
- handles different ploidy levels

APPLICATIONS:

- in screening studies
- to identify rare alleles and mutations
- to identify novel genetic markers
- to evaluate results obtained with other software's

